

INFUSION SET *specifications*


Right Angle, ClearView Right Angle Sets

RA Series

Right Angle Block Huber Sets

Right Angle Huber Infusion Sets feature the deflected/offset 'B' bevel point needle and a grip block bonded to a Tygon extension with a female luer lock and clamp. Grip blocks are color coded for easy needle gauge identification.

Right angle sets ship sterile (EtO) in boxes of 12. The RA series is available in a variety of gauges and needle and extension lengths. These sets are available with injection sites, microclave connectors, or caps.



CVRA Series

Right Angle Huber Disk Sets

ClearView Right Angle Infusion Sets feature the deflected/offset 'B' bevel point needle attached to a clear flexible disk with a female luer lock and clamp. The advantage of the disk is that it lies very flat against the skin improving needle stability. ClearView sets ship sterile (EtO) in boxes of 12. The CVRA series is available in a variety of gauges and needle and extension lengths. Sets are available with injection sites, microclave connectors, or caps.



Ordering Information

1-847-674-7131 (U.S.A.)

Needle Gauge	Tubing Length	Needle Length				
		3/8"	1/2"	5/8"	3/4"	1"
Right Angle Huber Infusion Sets - RA Series - with female luer and clamp						
22	6"/15cm	RA22-375-6	RA22-5-6	RA22-625-6	RA22-75-6	RA22-100-6
22	12"/30cm	RA22-375-12	RA22-5-12	RA22-625-12	RA22-75-12	RA22-100-12
22	18"/45cm	RA22-375-18	RA22-5-18	RA22-625-18	RA22-75-18	RA22-100-18
20	6"/15cm	RA20-375-6	RA20-5-6	RA20-625-6	RA20-75-6	RA20-100-6
20	12"/30cm	RA20-375-12	RA20-5-12	RA20-625-12	RA20-75-12	RA20-100-12
20	18"/45cm	RA20-375-18	RA20-5-18	RA20-625-18	RA20-75-18	RA20-100-18
19	6"/15cm	RA19-375-6	RA19-5-6	RA19-625-6	RA19-75-6	RA19-100-6
19	12"/30cm	RA19-375-12	RA19-5-12	RA19-625-12	RA19-75-12	RA19-100-12
19	18"/45cm	RA19-375-18	RA19-5-18	RA19-625-18	RA19-75-18	RA19-100-18
Right Angle Huber Disk Infusion Sets - CVRA Series - with female luer and clamp						
22	6"/15cm	CVRA22-375-6	CVRA22-5-6	CVRA22-625-6	CVRA22-75-6	CVRA22-100-6
22	12"/30cm	CVRA22-375-12	CVRA22-5-12	CVRA22-625-12	CVRA22-75-12	CVRA22-100-12
22	18"/45cm	CVRA22-375-18	CVRA22-5-18	CVRA22-625-18	CVRA22-75-18	CVRA22-100-18
20	6"/15cm	CVRA20-375-6	CVRA20-5-6	CVRA20-625-6	CVRA20-75-6	CVRA20-100-6
20	12"/30cm	CVRA20-375-12	CVRA20-5-12	CVRA20-625-12	CVRA20-75-12	CVRA20-100-12
20	18"/45cm	CVRA20-375-18	CVRA20-5-18	CVRA20-625-18	CVRA20-75-18	CVRA20-100-18
19	6"/15cm	CVRA19-375-6	CVRA19-5-6	CVRA19-625-6	CVRA19-75-6	CVRA19-100-6
19	12"/30cm	CVRA19-375-12	CVRA19-5-12	CVRA19-625-12	CVRA19-75-12	CVRA19-100-12
19	18"/45cm	CVRA19-375-18	CVRA19-5-18	CVRA19-625-18	CVRA19-75-18	CVRA19-100-18

Call for custom sizes